

**Plak-Vac is a reusable single patient oral hygiene instrument for patients with difficulty swallowing and expectorating, or who cannot perform routine oral care.**

The Plak-Vac Oral Suction Evacuator Brush provides the ultimate in care for patients requiring assistance with oral hygiene. A combination yankauer suction tube and toothbrush, the Plak-Vac allows the removal of debris, plaque, bacteria and fluid from the mouth of patients who are unable to swallow or have difficulty swallowing and expectorating. The Plak-Vac has soft bristles that are gentle to gums, teeth and other tissues in the mouth. Suction level is controlled by a finger control port located in the Plak-Vac handle.


Use the Plak-Vac to:

- Make oral hygiene simple and effective
- Remove mucus and food build-up, as well as plaque, bacteria and fluid
- Help eliminate chronic halitosis
- Reduce the risk of infection, aspiration and pneumonia
- Eliminate the need for separate yankauers, toothbrushes and suction catheters

Plak-Vac is ideally suited for use with both wall suction and portable pumps.

**Our Plak-Vac/Res-Q-Vac combination features the Plak-Vac oral suction evacuator brush with the Res-Q-Vac hand powered suction system.**

It is the perfect system for home care and long term care use where wall or portable suction may not be available, and since it does not require a power source it is completely portable.



Order #	Description
2200	Plak-Vac Oral Suction Evacuator Brush
2213	Plak-Vac/Res-Q-Vac Combination